SapphireAmp® Fast PCR Master Mix

Code No. RR350A
Size : 1 ml x 4
(for 160 PCR reactions)

Shipping at — 20°C
Store at — 20°C

Components : (for 160 reactions)
SapphireAmp® Fast PCR Master Mix (2X Premix) 1 ml×4
dH₂O 1 ml×4

Lot No.

Storage : — 20°C (at 4°C for 3 months)

Please store Master Mix at 4°C after thawing because the activity of the Master Mix may decrease by repeated freeze/thaws. Please mix gently and centrifuge briefly before each use.

Description :

SapphireAmp® Fast PCR Master Mix (2X Premix) is a premixed, 2-fold concentrated PCR reaction mixture suitable for amplification of DNA at the extension rate of 10 seconds/kb. This mixture is composed of a PCR enzyme for Hot Start PCR, dNTP Mixture, and buffer which are optimized for high-speed PCR. The master mix also contains a blue dye and a material for increasing the specific gravity of the solution thus enabling direct loading on an agarose gel for electrophoresis. PCR efficiency is greatly improved by the master mix format, simply add your primers, template, and begin the cycling process. Once the amplification is complete, the reactant can be applied directly to an agarose gel for electrophoresis.

This product is very useful for examining an insert fragment (recommended for less than 5 kb) by colony-directed PCR of *E. coli* recombinant cells, and for obtaining the result quickly (about 1 hour 15 minutes for 5 kb amplification). Furthermore, it is possible to amplify up to a 6 kb DNA fragment when using human genomic DNA as a template.

Application :

- DNA amplification by Polymerase Chain Reaction (PCR)
- Colony PCR

PCR test :

Good enzyme performance was confirmed by PCR amplification using human genome DNA (amplified fragment : 6 kb) as templates.

General reaction mixture for PCR (Total 50 µl) :

SapphireAmp® Fast PCR Master Mix (2X Premix)	25 µl
Forward Primer	0.2 µM (final conc.)
Reverse Primer	0.2 µM (final conc.)
Template	50 - 100 ng
plasmid DNA	100 pg - 10 ng
dH ₂ O (Sterilized distilled water)	up to 50 µl

PCR product :

As most PCR products amplified with SapphireAmp® Fast PCR Master Mix have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into a T-vector. It is also possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

Dye marker at electrophoresis :

When 5 µl of a reactant is loaded on a 1% TaKaRa Agarose L03 gel for electrophoresis, the blue markers are detected at positions corresponding to 1 kb and 3-5 kb. The dyes have absorptions at around 260nm and 620nm, respectively. The dyes may be removed by cutting out the gel or extracting DNA by NucleoSpin Extract II, if necessary.

PCR Condition (example) :

3 step PCR (human or mouse genomic DNA : <2 kb, bacterial genome DNA : <6 kb, colony PCR : <6 kb)

94°C	1 min.	→	98°C	5 sec.	} 30 cycles
			55°C	5 sec.	
			72°C	10 sec./kb *	
* : 5 sec. for less than 1 kb.					

2 step PCR (human or mouse genomic DNA 2 - 6 kb)

94°C	1 min.	→	98°C	5 sec.	} 30 cycles
			68°C	30 sec./kb	

Application :

A Mouse cDNA library, with insert lengths 0.5 - 5 kb, was constructed with pUC118. Colony-directed PCR was performed to check inserts.

- (1) The PCR reaction mixture was prepared as shown below, then 50 µl of the mixture was dispensed in each PCR tube.

For reaction of 8 samples :

	One reaction	For 9 reactions
SapphireAmp® Fast PCR Master Mix (2X Premix)	25 µl	225 µl
M13 Primer M4 (20 µM)	0.5 µl	4.5 µl
M13 Primer RV (20 µM)	0.5 µl	4.5 µl
dH ₂ O	24 µl	216 µl


- (2) A very small amount of the recombinant colony was picked with a sterilized micropipette tip, and was then suspended in the PCR reaction mixture that was prepared and dispensed into PCR tubes. The PCR reaction was performed following the condition below.

PCR condition :

94°C	1 min.	→	98°C	5 sec.	} 30 cycles (PCR : 1 hr 11 min.
			55°C	5 sec.	
			72°C	40 sec.	
					TaKaRa PCR Thermal Cycler Dice®)

Electrophoresis :

The PCR products were confirmed by electrophoresis.

1 - 8 : Colony PCR
M : 250 bp DNA Ladder (Dye Plus)
Load 5 µl of the reaction.
1% Agarose L03

NOTICE TO PURCHASER: LIMITED LICENSE

[L15] Hot Start PCR

Licensed under U.S. Patent No. 5,338,671 and 5,587,287 and corresponding patents in other countries.

[M57] LA Technology

This product is covered by the claims 6-16 of U.S. Patent No. 5,436,149 and its foreign counterpart patent claims.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com.

SapphireAmp® Fast PCR Master Mix

Code No. RR350A
Size : 1 ml×4
(for 160 PCR reactions)

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