




All outer dimensions include corner guards
Cage Stack Heights: Mouse-10; Large Mouse-7; Rat-6



72.5"W x 33"D x 82.3"H*



62.4"W x 33"D x 82.3"H*




54.2"W x 33"D x 82.3"H*


Type of Cage	Frame Part Number
Gentle-Air 160 cages	MAC-GA-F10D-16
Flex-Air 140 cages	MAC-FA-F10D-14

Type of Cage	Frame Part Number
Gentle-Air 140 cages	MAC-GA-F10D-14
Flex-Air 120 cages	MAC-FA-F10D-12


Type of Cage	Frame Part Number
Gentle-Air 120 cages	MAC-GA-F10D-12
Flex-Air 100 cages	MAC-FA-F10D-10



72.5"W x 26"D x 82.3"H*



62.4"W x 26"D x 82.3"H*




54.2"W x 26"D x 82.3"H*


Type of Cage	Frame Part Number
Gentle-Air 80 cages	MAC-GA-F10S-8
Flex-Air 70 cages	MAC-FA-F10S-7

Type of Cage	Frame Part Number
Gentle-Air 70 cages	MAC-GA-F10S-7
Flex-Air 60 cages	MAC-FA-F10S-6

Type of Cage	Frame Part Number
Gentle-Air 60 cages	MAC-GA-F10S-6
Flex-Air 50 cages	MAC-FA-F10S-5




72.5"W x 29"D x 82.3"H*



62.4"W x 29"D x 82.3"H*

Type of Cage	Frame Part Number
Gentle-Air 80 cages	MAC-GA-F10S-8-MR
Flex-Air 70 cages	MAC-FA-F10S-7-MR
Rat 36 cages	MAC-RV-F6S-6-MR
Lg. Mouse 42 cages	MAC-RV-F6S-6-MR

Type of Cage	Frame Part Number
Gentle-Air 70 cages	MAC-GA-F10S-7-MR
Flex-Air 60 cages	MAC-FA-F10S-6-MR
Rat 30 cages	MAC-FA-F6S-5-MR
Lg. Mouse 35 cages	MAC-RV-F6S-6-MR





GA - 39.29"W x 39.29"D x 82.75"H**

FA - 42.75"W x 42.75"D x 82.75"H**

Type of Cage	Frame Part Number
Gentle-Air 120 cages	MAC-GA-F10S-12 58" Diameter
Flex-Air 120 cages	MAC-FA-F10S-12 61" Diameter

*Frame Dimensions Without AHU. With AHU - 94.0"H
**Frame Dimensions Without AHU. With AHU - 96.2"H

Wall-Mount Frames / Cage Stack Heights: Mouse -10; Large Mouse-7; Rat-6

Main Frames that Will Support an AHU	Frame Dimensions	Type of Cage	Number of Cages	Frame Part Number
	67.125"W x 19"D x 93.0"H	Gentle-Air	80 cages	MAC-GA-F10W-8
	67.125"W x 19"D x 93.0"H	Flex-Air	70 cages	MAC-FA-F10W-7
	67.125"W x 23.5"D x 93.0"H	Rat Lg. Mouse	30 cages 35 cages	MAC-RV-F6W-5-MR
	59"W x 19"D x 93.0"H	Gentle-Air	70 cages	MAC-GA-F10W-7
	59"W x 19"D x 93.0"H	Flex-Air	60 cages	MAC-FA-F10W-6
	59"W x 23.5"D x 93.0"H	Rat Lg. Mouse	30 cages 35 cages	MAC-RV-F6W-5-MRS
	50.875"W x 19"D x 93.0"H	Gentle-Air	60 cages	MAC-GA-F10W-6
	50.875"W x 19"D x 93.0"H	Flex-Air	50 cages	MAC-FA-F10W-5
	50.875"W x 23.5"D x 93.0"H	Rat Lg. Mouse	24 cages 28 cages	MAC-RV-F6W-4-MR
Expansion Frames that Connect to Main Frames for Air Supply	Frame Dimensions	Type of Cage	Number of Cages	Frame Part Number
	42.75"W x 19"D x 93.0"H	Gentle-Air	50 cages	MAC-GA-F10W-5
	42.75"W x 19"D x 93.0"H	Flex-Air	40 cages	MAC-FA-F10W-4
	42.75"W x 23.5"D x 93.0"H	Rat Lg. Mouse	18 cages 21 cages	MAC-RV-F6W-3-MR
	34.625"W x 19"D x 93.0"H	Gentle-Air	40 cages	MAC-GA-F10W-4
	34.625"W x 19"D x 93.0"H	Flex-Air	30 cages	MAC-FA-F10W-3
	34.625"W x 23.5"D x 93.0"H	Rat Lg. Mouse	12 cages 14 cages	MAC-RV-F6W-2-MR
	25.6"W x 19"D x 93.0"H	Gentle-Air	30 cages	MAC-GA-F10W-3
	25.6"W x 19"D x 93.0"H	Flex-Air	20 cages	MAC-FA-F10W-2
	25.6"W x 23.5"D x 93.0"H	Rat Lg. Mouse	12 cages 14 cages	MAC-RV-F6W-2-MRS