



ALTERNATIVE DESIGN
MANUFACTURING & SUPPLY, INC.

Proven Products. Creative Solutions.

Material Handling & Storage Products



Accessories

Accessory Basket	3
Low Profile Water Bottle Basket	3
Shelf Rack Adjustable - Extra Shelf	4
Universal Water Bottle Basket	4
Water Bottle Basket	5
Water Bottle Basket 24 Position	5

Carts

Autoclave Cart with Adjustable Shelves	6
Bedding Cart	6
Bulk Truck Adjustable Shelf	7
Bulk Truck Cover	7
Bulk Truck Grid Doors - Half/Full	8
Bulk Truck Heavy Duty	8
Bulk Truck Medium Duty	9
Bulk Truck - Mini	9
Pan Wash & Storage Cart	10
Platform Cart - 2 Handle	10
Utility Cart - 2 Shelf	11
Water Bottle Basket Dollie	11

Racks

Cage Wash Rack	12
Cage Wash Rack - Open Design	12
Dunnage Rack	13
Shelf Rack -Adjustable	13
Shelf Rack - Fixed	14
Shelf Rack with Wire Shelves - Adjustable	14
Pan Wash Rack	15

Tables

Clean Table	15
-------------	----

For more information about the features and specs, including part numbers and product dimensions on each individual product, visit our website. alternativedesign.com

Overview


Features

Specs

Care

Warranty





Standard Features Include:

- Made from electro polished, 304 stainless steel
- Frame made from 1" angled steel for greater strength and durability
- 1/2" x 1/2" wire mesh sides and bottom
- 2 easy grab handles
- Autoclavable

Available Options:

- Accessory Basket only (no lid)
- Accessory Basket with lid
- Accessory Basket with latching lid (secure fit spring latch)

Part Numbers:

AB-18x12.25x6 (Basket Only)
AB-18x12.25x6-VL (Basket with Lid)
AB-18x12.25x6-L (Lid Only)

Dimensions:

18" W x 12.25" D x 8" H
18" W x 12.25" D x 8" H
18" W x 12.25" D x 8" H

Low Profile Water Bottle Basket



Standard Features Include:

- Lower height for better ergonomics
- 304 stainless steel wire frame
- Holds up to 20, 16 oz reduced height water bottles
- Autoclavable

Available Options:

- Low Profile Bottle Basket only (no lid)
- Low Profile Bottle Basket with lid
- Low Profile Bottle Basket with latching lid (secure fit spring latch)

Part Numbers:

WBB20-16RH (No Lid)
WBB20-16RH-VL (w/Lid)
WBB20-16RH-L (Latching Lid)

Dimensions:

18.5" W x 12.5" D x 8.5" H
18.5" W x 12.5" D x 8.5" H
18.5" W x 12.5" D x 8.5" H



Standard Features Include:

- Add an extra shelf to your adjustable shelf rack
- Heavy duty construction for durability and long life
- Made from 304 stainless steel
- Autoclavable

Part Numbers:

SRASH60X24
SRASH60X24ASM (w/Assembly)

Dimensions:

60" W x 24" D
60" W x 24" D

Universal Water Bottle Basket



Standard Features Include:

- Lower height for better ergonomics
- 304 stainless steel wire frame
- Holds up to 20, 16 oz reduced height water bottles
- Autoclavable

Available Options:

- Universal Bottle Basket only (no lid)
- Universal Bottle Basket with lid
- Universal Bottle Basket with latching lid (secure fit spring latch)

Part Numbers:

WBB24-U (No Lid)
WBB24-U-WL (w/Lid)
WBB24-U-L (Latching Lid)

Dimensions:

18.5" W x 12.5" D x 9" H
18.5" W x 12.5" D x 9" H
19" W x 12.5" D x 9" H



Standard Features Include:

- 304 stainless steel wire frame
- Holds up to 20 reduced height water bottles
- Autoclavable

Available Options:

- Bottle Basket only (no lid)
- Bottle Basket with lid
- Bottle Basket with latching lid (secure fit spring latch)

Part Numbers:

WBB20-16 (No Lid)
WBB20-16-WL (w/Lid)
WBB20-16-L (Latching Lid)

Dimensions:

18.5" W x 12.5" D x 13.75" H
18.5" W x 12.5" D x 13.75" H
18.5" W x 12.5" D x 13.75" H

Water Bottle Basket 24 Position



Standard Features Include:

- Holds up to 24, 8 oz french square, or 9 oz reduced height bottles
- 304 stainless steel wire frame
- Autoclavable

Part

WBB24-9RH-WL

Dimensions:

18.5" W x 12.5" D x 8.5" H

**Standard Features Include:**

- 304 stainless steel construction
- Perforated base
- Includes 3 perforated shelves
- Shelves can be placed in 10 different positions
- Locking bar allows two carts to be joined together
- Swivel casters, 2 with step brakes
- Holds up to 500 lbs

Part #

AUTOCLAVE CART

Dimensions:

36 3/4" W x 24 3/4" D x
57 5/8" H

Bedding Cart

**Standard Features Include:**

- 304 stainless steel construction
- Adjustable height (from 36" to 42") with quick release pins
- Includes 2 stainless steel shelves and 45 gallon Rubbermaid® container (42" L x 21" W)
- Stack cages on one of 3 collapsible stainless steel side trays
- One large side tray (40" L x 18" W) and two medium side trays (19" L x 18" W)
- Four swivel casters, 2 with brakes

Part Numbers:

BC-42X21X36-45G-V2 (Three Side Trays)
BC-42X21X36-45G (Two Side Trays)

Dimensions:

42"W x 21"D x 36"H
42"W x 21"D x 36"H



Standard Features Include:

- Heavy duty construction for durability and long life
- 304 stainless steel
- Autoclavable

Part Numbers:

BTS-60X24
BTS-60X32
BTS-44X27

Dimensions:

60"W x 24"D
44"W x 32"D
44"W x 27"D



Standard Features Include:

- Fits 60" x 24" x 70" size bulk trucks
- Autoclavable
- Available with either reinforced zippers or hook & loop fasteners
- Custom sizes also available (please inquire)

Part Numbers:

EC-60X24X68-Z
EC-60X24X68-V
EC-60X32X68-Z
EC-60X32X68-V

Dimensions:

Fits size 60" x 24" x 70"
Fits size 60" x 24" x 70"
Fits size 60" x 32" x 70"
Fits size 60" x 32" x 70"

**Standard Features Include:**

- Heavy duty construction for durability and long life
- 304 stainless steel
- Autoclavable

Part Numbers:**Full**

BTD-60X68-F
BTD-44X75-F

Half

BTD-60X34-H
BTD-44X37-H

Dimensions:**Full**

60"W x 60"H
44"W x 67"H

Half

60"W x 33"H
44"W x 37"H

**Standard Features Include:**

- Rated to hold up to 1,500 lbs
- Heavy duty construction for durability and long life
- 304 stainless steel, 1.25" square tube frame with 3/8" rods
- Solid metal floor, sloped for water drainage
- Smooth surfaces for easy sterilization
- Large cage-in stacking surface
- 6" Stainless steel casters for easier transporting, 2 with brakes
- Autoclavable

Part Numbers:

BT-60X24X68
BT-60X32X68

Dimensions:

60"W x 24"D x 68"H
60"W x 32"D x 68"H

**Standard Features Include:**

- Rated to hold up to 1000 lbs
- All 304 stainless steel construction
- Square tube/C channel frame for strength
- Solid bottom shelf with zero standing water design
- 5" grease zerk casters, 2 with brakes
- Autoclavable

Part Numbers:

BT-CH-44X27X68-MD
 BT-CH-45X29X68-MD
 BT-CH-60X24X68-MD
 BT-CH-60X32X68-MD

Dimensions:

44"W x 27"D x 68"H
 45"W x 29"D x 68"H
 60"W x 24"D x 68"H
 60"W x 32"D x 68"H

**Standard Features Include:**

- All 304 stainless steel construction
- Square tube/C channel frame for strength
- Grease zerk casters, 2 with brakes
- Autoclavable

Available Options:

- Removable shelf
- Full length door

Part #

BT-CH-36X24X36-MD

Dimensions:

36"W x 24"D x 36"H



Standard Features Include:

- Holds up to 22 metabolism pans
- 304 stainless steel construction
- Autoclavable
- 5" stainless steel swivel casters, 2 with brakes

Part

WRP-97X32X53-C

Dimensions:

97" W x 32" D x 53.4" H

Platform Cart - 2 Handle



Standard Features Include:

- 304 stainless steel construction
- Handle bars at both ends for convenience
- Center and caster braces for greater support
- Bed is 14 gauge for durability
- 6" swivel stainless steel casters to move easily, 2 with brakes
- Holds up to 400 lbs

Part

PFC-43X24X35

Dimensions:

43" W x 24" D x 34.4" H



Standard Features Include:

- 304 stainless steel construction
- 1" tube uprights for greater strength
- Shelves are made from 18 gauge with a 2B finish
- 5" swivel stainless steel casters for easy handling, 2 with brakes
- Holds up to 150 lbs

Available Options:

- 1 or 2 handlebars
- 3rd shelf
- 3 raised edges on the top shelf with 1 handlebar

Part Numbers:

UC2-40X24X36
 UC2-32X19X36
 UC2-40X24X36-1HB (Handlebars)
 UC2-32X19X36-1HB
 UC3-40X24X36 (3rd Shelf)
 UC2-40X24X36-1HB-3ES (Raised edges)

Dimensions:

40" W x 24" D x 36" H
 32" W x 19" D x 36" H
 40" W x 24" D x 36" H
 32" W x 19" D x 36" H
 40" W x 24" D x 36" H
 40" W x 24" D x 36" H

Water Bottle Basket Dollie



Standard Features Include:

- 304 stainless steel frame
- Off-set handle for easy handling
- 4" swivel stainless steel casters to move easily
- Solid welded frame
- Holds up to eight 9-inch high bottle baskets

Part

WBBD-6-SS

Dimensions:

20" W x 26.125" D x 36" H

**Standard Features Include:**

- Holds standard mouse, large mouse, and rat plastic cages
- 8 baskets made from 1" stainless mesh. Each basket will hold up to 10 mouse cages, 4 large mouse cages, or 4 rat cages
- 304 stainless steel construction
- 5" stainless steel swivel casters with grease zerks, 2 casters with brakes

Part #

WRC-76X29X70

Dimensions:

76" W x 29" D x 70" H

Cage Wash Rack - Open Design

**Standard Features Include:**

- Open design for more thorough cleaning
- Shelves auto adjust for all cage sizes. Spring loaded to assure tight fit
- 304 stainless steel construction
- 5" Stainless steel swivel casters, 2 with brakes

Part #

WRC-76X29X68-O

Dimensions:

76.1" W x 29" D x 68.5" H

**Standard Features Include:**

- Store items such as boxes and bags of feed, or bedding off the floor
- Made from 2" square tube extruded aluminum
- Rated to hold up to 4,000 lbs

Part #

DR-48X24X12-A

Dimensions:

48" W x 24" D x 12" H

Shelf Rack - Adjustable

**Standard Features Include:**

- Vertical uprights are made of 1" square tubing for greater strength from 304 stainless steel (#3 Finish, 16 Gauge)
- Shelves are made from 304 stainless steel (2B Finish, #18 Gauge)
- Swivel casters are 5" diameter by 1.25" tread *Maxim* wheel with expandable stem mount, grease zerks, and pedal wheel brakes on 2 of the casters
- Assembly adjusting bracket including a bolt plate, socket head screw, and washers
- Plastic end caps in upright
- Includes Allen wrenches needed for assembly
- Autoclavable

Part Numbers:

SRA7-60X24X72
 SRA6-60X24X72
 SRA5-60X24X72
 SRA4-60X24X72

Shelves:

7
 6
 5
 4

**Standard Features Include:**

- Vertical uprights are made of 1" square tubing for greater strength from 304 stainless steel (#3 Finish, 16 Gauge)
- Shelves are made from 304 stainless steel (2B Finish, #18 Gauge)
- Swivel casters are 5" diameter by 1.25" tread *Maxim* wheel with expandable stem mount, grease zerks, and pedal wheel brakes on 2 of the casters with stainless steel hardware
- Weight rated to hold up to 150 lbs. per shelf
- Autoclavable

Part Numbers:

SRF6-60X24X72
 SRF5-60X24X72
 SRF4-60X24X72

Shelves:

6
 5
 4

**Standard Features Include:**

- Rack has vertical uprights of 1" round stainless steel tubing with a 16-gauge wall
- Wire mesh shelves are adjustable in height and locked into place by setscrews
- Access to the set screws may be made from the outside edge of each shelf
- Rack is equipped with 5" standard expandable casters
- Two of the casters are equipped with step brakes
- Shipped unassembled and includes detailed assembly instructions
- Autoclavable

Part Numbers:

SRA5-60X24X72-VWS
 SRA5-60X24X66-VWS

Dimensions:

60"W x 24"D x 72"H
 60"W x 24"D x 66"H

**Standard Features Include:**

- Suitable for all types of excreta pans
- Includes adjustable pan dividers
- 304 stainless steel construction
- 5" stainless steel swivel casters, 2 with brakes

Part #

WRP-66X40X54

Dimensions:

66" W x 40" D x 54" H

**Standard Features Include:**

- 304 stainless steel construction with #3 finish
- 1" square tube legs for durability
- Leveler feet for stability
- 1-1/2" table top edge for greater rigidity
- Table surface with backsplash option is seamless for easier cleaning and less likely bacteria growth

Available Options:

- 1 or 2 handlebars
- 3rd shelf
- 3 raised edges on the top shelf with 1 handlebar

Part Numbers:

CT-60X24X34 (Standard)
 CT-60X24X34-B (With backsplash)
 CT-60X24X34-B-FS (With backsplash and full shelf)
 CT-48X24X34-FS (With backsplash and full shelf)

Dimensions:

60" W x 24" D x 34" H
 60" W x 24" D x 34" H
 60" W x 24" D x 34" H
 48" W x 24" D x 34" H

Alternative Design Manufacturing & Supply

Quality products and service for all your poultry and avian housing needs.

Whether it's an animal housing design to maximize your space usage, or material handling equipment to improve workflow efficiency, our commitment is to work closely with you to deliver the animal room solution that best meets your facility's space, needs, and budget.

Warranty



Alternative Design Manufacturing & Supply offers a 5-year warranty on all stainless steel fabricated products. These products are warrantied to be free from material and workmanship defects during this period. The warranty period will begin from either the date put into service, installation (if applicable), or 6 months from date of shipment, whichever occurs first.

- Individually Ventilated Caging (IVC) systems including rack frames and fabricated components
- All other stainless steel fabricated caging systems
- Fabricated material handling, storage and processing equipment made of stainless steel



Alternative Design offers a 2-year warranty on all Polysulfone plastic caging products and electrical components. These products are warrantied to be free from material and workmanship defects during this period. The warranty period will begin from the date of delivery. Our Polysulfone products will have a 2-year prorated warranty*. Please note – Polycarbonate caging is not recommended for autoclaving and any warranty will be voided if autoclaved.

- All electrical components
- Polysulfone plastic caging products or components

* Plastic caging materials that show physical, chemical, or heat damage will not be covered under warranty.



Alternative Design offers a 1-year warranty on all other products. These products are warrantied to be free from material and workmanship defects during this period. The warranty period will begin from the date of delivery unless otherwise stated.

- All other products manufactured or sold by Alternative Design Manufacturing & Supply

Should a failure occur within any of the the above listed time frames, Alternative Design will have the option to repair or replace the defective product(s), part(s), or component(s). This warranty is extended only to the original purchaser of the product and excludes damage due to abuse, misuse, as well as normal wear and tear. All non-warranty repairs will include a 90-day warranty from the date of delivery.

For more information on our industry-leading warranties, please contact Customer Service.



Warranty Repair or Service

It is the goal of Alternative Design to provide a quality product with full service warranty support. Ask about Alternative Design's Performance Maintenance Agreement program. You may obtain warranty repair or service by contacting our Customer Service Department.

Alternative Design Manufacturing & Supply, Inc.
3055 Cheri Whitlock Drive, P.O. Box 6330
Siloam Springs, AR 72761
USA

800-320-2459
479-524-4343
479-524-4125 (FAX)

info@altdesign.com
www.altdesign.com
f t in g+ YouTube