

MPC-5V312

2~8°C Pharmacy Refrigerator



MPC-5V312 2~8°C Pharmacy Refrigerator (312 liter) is widely used in hospitals pharmacy stores and laboratories. This is for placing and preserving storage such as drugs, vaccines, enzymes, hormones, stem cells, extracted RNA and gene libraries and important biological and chemical reagents.

FEATURERS



Fast pull down

Forced air cooling system



Heating glass

Double layer Low-E glass for fast defogging



Suction circulation system

Big fans for better temperature uniformity



Universal castor wheels

Proven heavy-duty lockable castors



Lockable doors

Standard lock with keys



High-quality materials

Antibacterial and corrosion resistant



Precise temp. control

Temp. control accurate to 0.1°C



5 types of alarm

Audible and visual alarms




Variable space

Height adjustable shelves




Eco-friendly



CFC-free refrigerant and energy-saving



MPC-5V312 Pull Down Curve



STRUCTURE



PARAMETER

Performance		Refrigeration System		Alarms	
Cooling method	Forced Air Cooling	Refrigerant	R290	High/low temperature	Y
Defrost	Auto	Temp. controller	Microprocessor	Power failure	Y
Shelf/drawer	Shelves/5	Sensor	NTC	Sensor error	Y
Temperature range	2°C~8°C	Display	Digital display	Controller failure	Y
Ambient temperature	10~32°C	External material	SPCC sprayed steel	Door ajar	Y
Noise (db)	50	Compressor QTY	1	Power failure backup	8 hours

Electrical Data		Dimensions		Other	
Voltage frequency(V/Hz)	220V 50Hz	Package size(W*D*H)(mm)	710*680*1960	Inside door	5 shelves
Power (W)	245	Capacity (L/cu.ft)	312/11	Remote alarm port	Optional
Current (A)	1.3	Exterior size(W*D*H)(mm)	625*594*1707	Access port QTY/Dia.	1/25mm
Power consumption	1.85kWh/24h	Interior size(W*D*H)(mm)	560*485*1350	Internal material	HIPS