

**Applied Biological Materials Inc.**

Telephone: 1-866-757-2414

Email: info@abmgood.comWebsite: www.abmGood.com

3M Sodium Acetate Solution, pH 5.2

Cat.No.	Unit
CH126	500 ml

Cat. No.	CH126
Name	3M Sodium Acetate Solution, pH 5.2
Unit	500 ml
Category	Buffers & Chemicals
Description	Sodium Acetate, a solution used in DNA extraction, reverse crosslinking and DNA purification, and to precipitate Gata6 (GATA binding protein 6) cDNA as a part of high-purity plasmid preparation.
Concentration	3M
Storage Condition	Room temperature

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information (1-866-757-2414).